

GW3-TRBO 2.15.5 Release Notes

New Features

- GW3-TRBO is now compatible with PCR 2.9 firmware
- IPSC Digital Backhaul repeaters are now displayed in GW3-TRBO
- Capacity Max telephone calls are now displayed and archived in GW3-TRBO
- Telephone call filtering options are now available in several reports
- Site and talkgroup affiliation information are now available in the GW3-TRBO REST API
- GW3-TRBO client software can be installed and updated via REST API connection to the Host
- GW3-TRBO clients can automatically download and install updates from the host when the host is upgraded
- GW3-TRBO clients can now be installed and upgraded via Active Directory Group Policy
- GEnSAC now unmutes console during console PTT allowing dispatcher to hear the talk permit tone or busy tone from the control station

Improvements and Fixes

- Improved logic for tabbing between fields in all products
- GW3-TRBO Alerts GUI notification window now remembers its last position and size
- GW3-TRBO and iVista products now use .NET 4.6.2
- GW3-TRBO busy reports now show in the appropriate report categories
- SQL Server no longer needs to be deployed in 'mixed mode' in order to install GW3-TRBO
- Busy Summary by Site report wizard now specifies a unit of time for the 'Max Duration' parameter
- GW3-TRBO now supports 100+ client connections
- Mass Random Access diagnostic messages are now displayed in GW3-TRBO and iVista
- Removed obsolete documentation and features from GenFIT product
- Alias will automatically remove the 'Offline' when it sees traffic on a channel
- iVista Map help documentation is improved
- iVista now supports 100+ client connections
- iVista Login Audit page overlapping text issue has been corrected
- iVista Users will be redirected to the Login Page upon session timeout
- iVista provisioning for Connect Plus no longer displays duplicate sync messages
- iVista provisioning for Connect Plus will now allow for duplicated radio aliases